

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

3 plane localiser - 3 plane localiser

IMAGING PARAMETERS	
Imaging Mode	2D
Pulse Sequence	Spin Echo
Imaging Options	Seq, EDR, Fast, SS, ARC
Phase	2.50
SCAN RANGE	
FOV	50.0
Slice Thickness	6.00
Slice Spacing	5.0
Overlap Locations	0
ACQ TIMING	
Freq	256
Phase	128
Freq DIR	Unswap
# of Acq. Before Pause	0
Phase FOV	1.00
Auto Shim	On
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective
FMRI	
PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

SCAN TIMING	
TE	100.0
Number of Echoes	1
TR	Minimum
Receiver Bandwidth	83.33
IMAGE ENHANCE	
Filter Choice	None
USER CVS	
User CV2	240.00
User CV13	1.00
User CV Mask2	256
MULTI-PHASE	
Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0
DIFFUSION	
Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;
CONTRAST	
Contrast Yes/No	No

3 plane localiser - 3 plane localiser

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

ARDL H Cor SSFSE 3mm TE 60 ARDLHi - ARDL H Cor SSFSE 3mm TE 60 ARDLHi	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	2D	TE	62.0
	Pulse Sequence	Spin Echo	Number of Echoes	1
	Imaging Options	EDR, TRF, Fast, SS, ARC	TR	Minimum
	Phase	2.00	Receiver Bandwidth	83.33
	SCAN RANGE		IMAGE ENHANCE	
	FOV	40.0	Filter Choice	None
	Slice Thickness	3.00	USER CVS	
	Slice Spacing	0.3	User CV1	1.00
	Overlap Locations	0	User CV Mask2	256
	Number of Slices	50	MULTI-PHASE	
	ACQ TIMING		Seperate Series	0
	Freq	288	Mask Phase	0
	Phase	224	Mask Pause	0
	Freq DIR	S/I	Preserve	0
	NEX	1.00	DIFFUSION	
	# of Acq. Before Pause	25	Recon All Images	On
	Phase FOV	1.00	# Synthetic b-values	1
	Auto Shim	Auto	Synthetic b-value	1000.0;
	Phase Correction	No	CONTRAST	
	RF Drive Mode	Single	Contrast Yes/No	No
	Excitation Mode	Selective		
ARDL H Cor SSFSE 3mm TE 60 ARDLHi - ARDL H Cor SSFSE 3mm TE 60 ARDLHi	FMRI			
	PSD Trigger	Internal		
	View Order	Bottom/Up		
	# of Repetitions REST	0		
	# of Repetitions ACTIVE	0		
	SAT			
	Tag Type	None		
	TRICKS			
	Pause On/Off	On		
	Auto Subtract	0		
	Auto SCIC	2		

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

ARDL H Ax Dual Echo ARDL ZIP 512 H - ARDL H Ax Dual Echo ARDL ZIP 512 H

IMAGING PARAMETERS

Imaging Mode 2D
Pulse Sequence SPGR
Imaging Options EDR, Fast, ZIP512, ARC
Phase 2.00

SCAN RANGE

FOV 40.0
Slice Thickness 5.00
Slice Spacing 0.6
Overlap Locations 0
Number of Slices 36

ACQ TIMING

Freq 224
Phase 224
Freq DIR R/L
NEX 1.00
of Acq. Before Pause 1
Phase FOV 0.80
Auto Shim Auto
Phase Correction No
RF Drive Mode Single
Excitation Mode Selective

USER CVS

User CV20 1.00
TR Min 120.0
TR Max 250.0
User CV Mask2 0

MULTI-PHASE

Seperate Series 0
Mask Phase 0
Mask Pause 0
Preserve 0

DIFFUSION

Number of Diffusion Directions 0
Dual Spin Echo Off
Diffusion Tensor No Selection
Processing Output
Recon All Images On
Synthetic b-values 1
Synthetic b-value 1000.0;

CONTRAST

Contrast Yes/No No

SCAN TIMING

Flip Angle 60
Number of Echoes 2
TR 120.0

IMAGE ENHANCE

Filter Choice None

GATING/TRIGGER

Pause After Navigator Prescan 0

FMRI

PSD Trigger Internal
View Order Bottom/Up
of Repetitions REST 0
of Repetitions ACTIVE 0

SAT

SAT Location S
SAT Location I
Tag Type None

TRICKS

Pause On/Off On
Auto Subtract 0
Auto SCIC 2

ARDL H Ax Dual Echo ARDL ZIP 512 H - ARDL H Ax Dual Echo ARDL ZIP 512 H

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

IMAGING PARAMETERS

Imaging Mode	3D
Pulse Sequence	LAVA
Imaging Options	EDR, Fast, ARC, HS
PSD Name	efgre3d_aspir
Phase	2.00
Slice	1.30
HyperSense	1.10

SCAN RANGE

FOV	44.0
Slice Thickness	2.20
Location per Slab	112
Overlap Locations	0
Number of Slices	1

ACQ TIMING

Freq	288
Phase	200
Freq DIR	R/L
NEX	1.00
Phase FOV	0.70
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
Fat/Water Saturation	Fat Special

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	2

SCAN TIMING

Flip Angle	12.0
Number of Echoes	1
TI	24
Receiver Bandwidth	50.00

IMAGE ENHANCE

Filter Choice	None
---------------	------

USER CVS

User CV4	10.20
User CV6	1.00
User CV34	90.00
User CV Mask2	8

MULTI-PHASE

Seperate Series	0
Trigger Delay without AV	10.2
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST


Contrast Yes/No	No
-----------------	----

Ax LAVA ASPIR ardl med pre - Ax LAVA ASPIR ardl med pre

Ax LAVA ASPIR ardl med pre - Ax LAVA ASPIR ardl med pre

OTHERS

Protocol Notes



LAVA Considerations:
-Use LAVA for abdominal scanning, in particular, liver imaging
-LAVA automatically turns on the ARC Imaging Option. If you deselect ARC and select ASSET, you must have acquired a calibration scan prior to acquiring the LAVA scan. If any or all of these scan parameters are selected (ASSET, Hyperband, PURE), and if you select ON from the Calibration in Prescan menu, which is located on the Details tab, a calibration scan is acquired during Auto Prescan
Consider the following with Fat/Water SAT:
-The overall SNR of Improved Fat SAT images may slightly decrease. Consider adjusting parameters with affect SNR, e.g. increase NEX, increase the FOV, decrease matrix size, decrease the Receive Bandwidth
-Uneven saturation can still occur as a result of local inhomogeneities, e.g. at air/tissue interfaces when the anatomy of interest is non-uniform
-Site or patient specific inhomogeneities may be unavoidable even at isocenter. The result can be uneven suppression.
Techniques for Fat Saturation include:
Chemical Fat SAT, STIR, ASPIR, and Flex.
-Consider using shom volumes to increase Fat SAT homogeneity.
Localized TG is used to improve the accuracy of the TG for a specifically designated area, such as cardiac, pelvis or prostate imaging, or when an anatomy is positioned at a patient's side. Localize TG can be used with any pulse sequence to improve signal homogeneity and is found under the Shim Volume menu

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

Ax LAVA ASPIR Mph ardl MED - Ax LAVA ASPIR Mph ardl MED

Ax LAVA ASPIR Mph ardl MED - Ax LAVA ASPIR Mph ardl MED

IMAGING PARAMETERS

Imaging Mode	3D
Pulse Sequence	LAVA
Imaging Options	EDR, Fast, MPh, SmartPrep, ARC, HS
PSD Name	efgre3d_aspir
Phase	2.00
Slice	1.30
HyperSense	1.10

SCAN RANGE

FOV	44.0
Slice Thickness	2.20
Location per Slab	112
Overlap Locations	0
Number of Slices	1

ACQ TIMING

Freq	288
Phase	200
Freq DIR	R/L
NEX	1.00
# of Acq. Before Pause	1
Phase FOV	0.70
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
Fat/Water Saturation	Fat Special

TRICKS

Pause On/Off	On
Auto Subtract	7
Auto SCIC	2

SCAN TIMING

Flip Angle	12.0
Number of Echoes	1
TI	24
Receiver Bandwidth	50.00

IMAGE ENHANCE

Filter Choice	None
---------------	------

USER CVS

User CV2	40.00
User CV4	4.00
User CV6	1.00
User CV34	90.00
User CV Mask2	8

MULTI-PHASE

Slice per Location	3
Delay after Acquisition	Minimum
Seperate Series	0
Delay after Acquisition without AV	17
Trigger Delay without AV	4
Mask Phase	0
Mask Pause	0

DIFFUSION


Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST

Contrast Yes/No	No
-----------------	----

OTHERS

Protocol Notes



LAVA Considerations:
-Use LAVA for abdominal scanning, in particular, liver imaging
-LAVA automatically turns on the ARC Imaging Option. If you deselect ARC and select ASSET, you must have acquired a calibration scan prior to acquiring the LAVA scan. If any or all of these scan parameters are selected (ASSET, Hyperband, PURE), and if you select ON from the Calibration in Prescan menu, which is located on the Details tab, a calibration scan is acquired during Auto Prescan
Consider the following with Fat/Water SAT:
-The overall SNR of Improved Fat SAT images may slightly decrease.
Consider adjusting parameters with affect SNR, e.g. increase NEX, increase the FOV, decrease matrix size, decrease the Receive Bandwidth
-Uneven saturation can still occur as a result of local inhomogeneities, e.g. at air/tissue interfaces when the anatomy of interest is non-uniform
-Site or patient specific inhomogeneities may be unavoidable even at isocenter. The result can be uneven suppression.
Techniques for Fat Saturation include:
Chemical Fat SAT, STIR, ASPIR, and Flex.
-Consider using shom volumes to increase Fat SAT homogeneity.
Localized TG is used to improve the accuracy of the TG for a specifically designated area, such as cardiac, pelvis or prostate imaging, or when an anatomy is positioned at a patient's side. Localize TG can be used with any pulse sequence to improve signal homogeneity and is found under the Shim Volume menu

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

ADRL Ax SSFSE TE 160 5mm ARDL med - ADRL Ax SSFSE TE 160 5mm ARDL med	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	2D	TE	162.0
	Pulse Sequence	Spin Echo	Number of Echoes	1
	Imaging Options	EDR, Fast, SS, ARC	TR	Minimum
	Phase	2.00	Receiver Bandwidth	90.91
	SCAN RANGE		IMAGE ENHANCE	
	FOV	42.0	Filter Choice	None
	Slice Thickness	5.00	USER CVS	
	Slice Spacing	0.3	User CV1	1.00
	Overlap Locations	0	User CV2	240.00
	Number of Slices	37	User CV13	1.00
	ACQ TIMING		User CV Mask2	256
	Freq	260	MULTI-PHASE	
	Phase	224	Seperate Series	0
	Freq DIR	R/L	Mask Phase	0
	NEX	1.00	Mask Pause	0
	# of Acq. Before Pause	16	Preserve	0
	Phase FOV	0.80	DIFFUSION	
	Auto Shim	Auto	Recon All Images	On
	Phase Correction	No	# Synthetic b-values	1
	RF Drive Mode	Single	Synthetic b-value	1000.0;
	Excitation Mode	Selective	CONTRAST	
ADRL Ax SSFSE TE 160 5mm ARDL med - ADRL Ax SSFSE TE 160 5mm ARDL med	FMRI		Contrast Yes/No	Yes
	PSD Trigger	Internal		
	View Order	Bottom/Up		
	# of Repetitions REST	0		
	# of Repetitions ACTIVE	0		
	SAT			
	Tag Type	None		
	TRICKS			
	Pause On/Off	On		
	Auto Subtract	0		
	Auto SCIC	2		

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

ARDL Ax SSFSE FS TE 90 3mm ARDL HIGH - ARDL Ax SSFSE FS TE 90 3mm ARDL HIGH

IMAGING PARAMETERS

Imaging Mode	2D
Pulse Sequence	Spin Echo
Imaging Options	EDR, Fast, SS, ARC
Phase	2.00

SCAN RANGE

FOV	42.0
Slice Thickness	5.00
Slice Spacing	0.5
Overlap Locations	0
Number of Slices	36

ACQ TIMING

Freq	260
Phase	256
Freq DIR	R/L
NEX	1.00
# of Acq. Before Pause	16
Phase FOV	0.80
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
Fat/Water Saturation	Fat Special

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	2

SCAN TIMING

TE	90.0
Number of Echoes	1
TR	Minimum
TI	130
Receiver Bandwidth	83.33

IMAGE ENHANCE

Filter Choice	None
---------------	------

USER CVS

User CV1	1.00
User CV2	240.00
User CV13	1.00
User CV Mask2	256

MULTI-PHASE

Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Number of Diffusion Directions	0
Dual Spin Echo	Off
Diffusion Tensor Processing Output	No Selection
Recon All Images	On
Multi b-values	1000.0;
Multi NEX Values	1.0;
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST

Contrast Yes/No	Yes
-----------------	-----

ARDL Ax SSFSE FS TE 90 3mm ARDL HIGH - ARDL Ax SSFSE FS TE 90 3mm ARDL HIGH

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

DWI DL Ax Nav high - DWI DL Ax Nav high	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	2D	TE	Minimum
	Pulse Sequence	Spin Echo	Number of Echoes	1
	Imaging Options	EDR, EPI, DIFF, Asset, Nav	Number of Shots	1
	Phase	2.00	IMAGE ENHANCE	
	SCAN RANGE		Filter Choice	None
	FOV	40.0	GATING/TRIGGER	
	Slice Thickness	5.00	Resp. Trigger Window	30
	Slice Spacing	0.4	Acceptance Window	2.0
	Overlap Locations	0	Run RR Measurement	1
	Number of Slices	37	Auto Trigger Time	5
	ACQ TIMING		Navigator Pre Scan Time	10
	Freq	90	Max RR Measurement Time	20
	Phase	100	FMRI	
	Freq DIR	R/L	PSD Trigger	Internal
	Phase FOV	1.00	View Order	Bottom/Up
	Auto Shim	Auto	# of Repetitions REST	0
	Phase Correction	Yes	# of Repetitions ACTIVE	0
	RF Drive Mode	Single	SAT	
	Excitation Mode	Selective	Tag Type	None
	USER CVS		TRICKS	
	User CV0	1.00	Pause On/Off	On
	User CV Mask2	256	Auto Subtract	0
	MULTI-PHASE		Auto SCIC	2
	Seperate Series	0	CONTRAST	
	Mask Phase	0	Contrast Yes/No	Yes
	Mask Pause	0		
	Preserve	0		
	DIFFUSION			
	Optimized TE	Yes		
	Diffusion Directions	All		
	Number of Diffusion Directions	3		
	Number of T2 Images	0		
	Dual Spin Echo	Off		
	Recon All Images	On		
	Multi b-values	100.0;600.0;		
	Multi NEX Values	2.0;4.0;		
	Real Time Field Adjustment	1		
	# Synthetic b-values	1		
	Synthetic b-value	800.0;		
	TRACKER			
	Tracker Length	120.0		
	Tracker Thickness	20.0		

DWI DL Ax Nav high - DWI DL Ax Nav high

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

3D Ax LAVA Star - 3D Ax LAVA Star

IMAGING PARAMETERS	
Imaging Mode	3D
Pulse Sequence	LAVA
Imaging Options	EDR, Fast, ARC, Star
PSD Name	efgre3d_aspir
Phase	1.50
Slice	2.00
SCAN RANGE	
FOV	40.0
Slice Thickness	2.20
Location per Slab	100
Overlap Locations	0
Number of Slices	1
ACQ TIMING	
Freq	288
Freq DIR	R/L
NEX	1.80
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective
FMRI	
PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0
SAT	
Tag Type	None
Fat/Water Saturation	Fat Special
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	2
OTHERS	
Protocol Notes	Streak artifacts can increase from arm movement and incomplete fat suppression To decrease streak artifacts: Increase NEX Decrease phase acceleration Raise arms above head Increase number of spokes=less streaks

SCAN TIMING	
Flip Angle	12
Number of Echoes	1
TI	23
Receiver Bandwidth	62.50
IMAGE ENHANCE	
Filter Choice	A
USER CVS	
User CV30	1.00
User CV34	70.00
User CV Mask2	8
MULTI-PHASE	
Seperate Series	0
Trigger Delay without AV	0
Mask Phase	0
Mask Pause	0
Preserve	0
DIFFUSION	
Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;
CONTRAST	
Contrast Yes/No	No

3D Ax LAVA Star - 3D Ax LAVA Star

3D Ax LAVA Star - 3D Ax LAVA Star

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

Ax IDEAL IQ - Ax IDEAL IQ	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	3D	Flip Angle	8
	Pulse Sequence	IDEAL IQ	TE	Min Full
	Imaging Options	EDR, Fast, IDEAL, ARC	Number of Echoes	6
	IDEAL	207	Echo Train Length	6
	Phase	2.00	Number of Shots	1
	Slice	1.00	Receiver Bandwidth	100.00
	SCAN RANGE		IMAGE ENHANCE	
	FOV	40.0	Filter Choice	None
	Slice Thickness	6.0	USER CVS	
	Location per Slab	32	User CV4	8.70
	Overlap Locations	0	User CV30	1.00
	Number of Slices	1	User CV Mask2	0
	ACQ TIMING		MULTI-PHASE	
	Freq	224	Seperate Series	0
	Phase	192	Trigger Delay without AV	8.7
	Freq DIR	R/L	Mask Phase	0
	NEX	0.50	Mask Pause	0
	Phase FOV	0.80	Preserve	0
	Auto Shim	Auto	DIFFUSION	
	Phase Correction	No	Recon All Images	On
	RF Drive Mode	Single	# Synthetic b-values	1
	Excitation Mode	Selective	Synthetic b-value	1000.0;
Ax IDEAL IQ - Ax IDEAL IQ	FMRI		CONTRAST	
	PSD Trigger	Internal	Contrast Yes/No	No
	Slice Order	Interleaved		
	View Order	Bottom/Up		
	# of Repetitions REST	0		
	# of Repetitions ACTIVE	0		
	SAT			
	Tag Type	None		
	TRICKS			
	Pause On/Off	On		
	Auto Subtract	0		
	Auto SCIC	Off		

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

Ax LAVA ASPIR ardl med delayed - Ax LAVA ASPIR ardl med delayed

IMAGING PARAMETERS

Imaging Mode 3D
Pulse Sequence LAVA
Imaging Options EDR, Fast, ARC, HS
PSD Name efgre3d_aspir
Phase 2.00
Slice 1.30
HyperSense 1.10

SCAN RANGE

FOV 44.0
Slice Thickness 2.20
Location per Slab 112
Overlap Locations 0
Number of Slices 1

ACQ TIMING

Freq 288
Phase 200
Freq DIR R/L
NEX 1.00
Phase FOV 0.70
Auto Shim Auto
Phase Correction No
RF Drive Mode Single
Excitation Mode Selective

FMRI

PSD Trigger Internal
View Order Bottom/Up
of Repetitions REST 0
of Repetitions ACTIVE 0

SAT

Tag Type None
Fat/Water Saturation Fat Special

TRICKS

Pause On/Off On
Auto Subtract 0
Auto SCIC 2

SCAN TIMING

Flip Angle 12.0
Number of Echoes 1
TI 24
Receiver Bandwidth 50.00

IMAGE ENHANCE

Filter Choice None

USER CVS

User CV4 10.20
User CV6 1.00
User CV34 90.00
User CV Mask2 8

MULTI-PHASE

Seperate Series 0
Trigger Delay without AV 10.2
Mask Phase 0
Mask Pause 0
Preserve 0

DIFFUSION

Recon All Images On
Synthetic b-values 1
Synthetic b-value 1000.0;


CONTRAST

Contrast Yes/No No

Ax LAVA ASPIR ardl med delayed - Ax LAVA ASPIR ardl med delayed

OTHERS

Protocol Notes



LAVA Considerations:
-Use LAVA for abdominal scanning, in particular, liver imaging
-LAVA automatically turns on the ARC Imaging Option. If you deselect ARC and select ASSET, you must have acquired a calibration scan prior to acquiring the LAVA scan. If any or all of these scan parameters are selected (ASSET, Hyperband, PURE), and if you select ON from the Calibration in Prescan menu, which is located on the Details tab, a calibration scan is acquired during Auto Prescan
Consider the following with Fat/Water SAT:
-The overall SNR of Improved Fat SAT images may slightly decrease.
Consider adjusting parameters with affect SNR, e.g. increase NEX, increase the FOV, decrease matrix size, decrease the Receive Bandwidth
-Uneven saturation can still occur as a result of local inhomogeneities, e.g. at air/tissue interfaces when the anatomy of interest is non-uniform
-Site or patient specific inhomogeneities may be unavoidable even at isocenter. The result can be uneven suppression.
Techniques for Fat Saturation include:
Chemical Fat SAT, STIR, ASPIR, and Flex.
-Consider using shom volumes to increase Fat SAT homogeneity.
Localized TG is used to improve the accuracy of the TG for a specifically designated area, such as cardiac, pelvis or prostate imaging, or when an anatomy is positioned at a patient's side. Localize TG can be used with any pulse sequence to improve signal homogeneity and is found under the Shim Volume menu

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

COR LAVA ASPIR ardl med - COR LAVA ASPIR ardl med

IMAGING PARAMETERS

Imaging Mode	3D
Pulse Sequence	LAVA
Imaging Options	EDR, Fast, ARC, HS
PSD Name	efgre3d_aspir
Phase	2.00
Slice	1.40
HyperSense	1.20

SCAN RANGE

FOV	45.0
Slice Thickness	2.30
Location per Slab	112
Overlap Locations	0

ACQ TIMING

Freq	300
Phase	192
Freq DIR	S/I
NEX	1.00
Phase FOV	0.80
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
Fat/Water Saturation	Fat Special

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	2

SCAN TIMING

Flip Angle	12.0
Number of Echoes	1
TI	24
Receiver Bandwidth	62.50

IMAGE ENHANCE

Filter Choice	None
---------------	------

USER CVS

User CV4	10.20
User CV6	1.00
User CV34	90.00
User CV Mask2	8

MULTI-PHASE

Seperate Series	0
Trigger Delay without AV	10.2
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Recon All Images	On
Multi b-values	1000.0;
Multi NEX Values	1.0;
# Synthetic b-values	1
Synthetic b-value	1000.0;


CONTRAST

Contrast Yes/No	No
-----------------	----

COR LAVA ASPIR ardl med - COR LAVA ASPIR ardl med

OTHERS

Protocol Notes



LAVA Considerations:
-Use LAVA for abdominal scanning, in particular, liver imaging
-LAVA automatically turns on the ARC Imaging Option. If you deselect ARC and select ASSET, you must have acquired a calibration scan prior to acquiring the LAVA scan. If any or all of these scan parameters are selected (ASSET, Hyperband, PURE), and if you select ON from the Calibration in Prescan menu, which is located on the Details tab, a calibration scan is acquired during Auto Prescan
Consider the following with Fat/Water SAT:
-The overall SNR of Improved Fat SAT images may slightly decrease.
Consider adjusting parameters with affect SNR, e.g. increase NEX, increase the FOV, decrease matrix size, decrease the Receive Bandwidth
-Uneven saturation can still occur as a result of local inhomogeneities, e.g. at air/tissue interfaces when the anatomy of interest is non-uniform
-Site or patient specific inhomogeneities may be unavoidable even at isocenter. The result can be uneven suppression.
Techniques for Fat Saturation include:
Chemical Fat SAT, STIR, ASPIR, and Flex.
-Consider using shom volumes to increase Fat SAT homogeneity.
Localized TG is used to improve the accuracy of the TG for a specifically designated area, such as cardiac, pelvis or prostate imaging, or when an anatomy is positioned at a patient's side. Localize TG can be used with any pulse sequence to improve signal homogeneity and is found under the Shim Volume menu

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

IMAGING PARAMETERS

Imaging Mode	2D
Pulse Sequence	Spin Echo
Imaging Options	Seq, EDR, TRF, Fast, SS, FR, ARC
Phase	2.00

SCAN RANGE

FOV	48.0
Slice Thickness	8.0
Slice Spacing	10.0

ACQ TIMING

Freq	384
Phase	160
Freq DIR	Unswap
# of Acq. Before Pause	0
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Quadrature
Excitation Mode	Selective

USER CVS

User CV1	1.00
User CV39	1.00
User CV Mask2	256

MULTI-PHASE

Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Recon All Images	On
Multi b-values	1000.0;
Multi NEX Values	1;
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST

Contrast Yes/No	No
-----------------	----

SCAN TIMING

TE	80.0
Number of Echoes	1
TR	700.0
Receiver Bandwidth	83.33

IMAGE ENHANCE

Filter Choice	None
---------------	------

GATING/TRIGGER

Auto Trigger Type	Off
-------------------	-----

FMRI

PSD Trigger	Internal
Slice Order	Interleaved
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
----------	------

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

IMAGING PARAMETERS

Imaging Mode *2D*
Pulse Sequence *Spin Echo*
Imaging Options *EDR, Fast, SS, ARC*
Phase *2.00*

SCAN RANGE

FOV *50.0*
Slice Thickness *4.00*
Slice Spacing *1.0*
Overlap Locations *0*
Number of Slices *20*

ACQ TIMING

Freq *256*
Phase *256*
Freq DIR *S/I*
NEX *1.00*
of Acq. Before Pause *15*
Phase FOV *1.00*
Auto Shim *Auto*
Phase Correction *No*
RF Drive Mode *Single*
Excitation Mode *Selective*

FMRI

PSD Trigger *Internal*
View Order *Bottom/Up*
of Repetitions REST *0*
of Repetitions ACTIVE *0*

SAT

Tag Type *None*

TRICKS

Pause On/Off *On*
Auto Subtract *0*
Auto SCIC *2*

SCAN TIMING

TE *90.0*
Number of Echoes *1*
TR *1000.0*
Receiver Bandwidth *90.91*

IMAGE ENHANCE

Filter Choice *None*

USER CVS

User CV1 *1.00*
User CV2 *240.00*
User CV13 *1.00*
User CV Mask2 *256*

MULTI-PHASE

Seperate Series *0*
Mask Phase *0*
Mask Pause *0*
Preserve *0*

DIFFUSION

Recon All Images *On*
Synthetic b-values *1*
Synthetic b-value *1000.0;*

CONTRAST

Contrast Yes/No *No*

sag SSFSE BH - sag SSFSE BH


sag SSFSE BH - sag SSFSE BH

OTHERS

Protocol Notes

General Considerations:

-To avoid contrast signal differences between contiguous slices, use a slice spacing $>20\%$ of slice thickness. If number of acquisitions is >1 , then 0 slice spacing can be used

*SSFSE Considerations:**-Use SSFSE:*

-to reduce motion artifact and imaging time
 -to scan uncooperative patients in short scan times
 -for breath hold abdominal and cardiac imaging
 -with long TE values (300-1300 ms) to image the gallbladder and biliary tree

*No Phase Wrap**Considerations:*

-When No Phase Wrap factor is set to a value greater than 1.0, it allows the patient to have arms at the side, which may increase patient comfort
 -A larger NPW value increases the ETL, which may increase image blurring. To preserve image sharpness decrease the value
 -A slight increase in SNR may occur due to more noise averaging that occurs with NPW

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

IMAGING PARAMETERS

Imaging Mode	3D
Pulse Sequence	LAVA
Imaging Options	EDR, Fast, ARC, HS
PSD Name	efgre3d_aspir
Phase	2.00
Slice	1.25
HyperSense	1.20

SCAN RANGE

FOV	42.0
Slice Thickness	1.30
Location per Slab	112
Overlap Locations	0
Number of Slices	1

ACQ TIMING

Freq	300
Phase	200
Freq DIR	S/I
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
Fat/Water Saturation	Fat Special

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	2

SCAN TIMING

Flip Angle	12.0
Number of Echoes	1
TI	24
Receiver Bandwidth	62.50

IMAGE ENHANCE

Filter Choice	None
---------------	------

USER CVS

User CV6	1.00
User CV34	90.00
User CV Mask2	8

MULTI-PHASE

Seperate Series	0
Trigger Delay without AV	0
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Recon All Images	On
Multi b-values	1000.0;
Multi NEX Values	1.0;
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST


Contrast Yes/No	No
-----------------	----

3D cor LAVA ASPIR BH - 3D cor LAVA ASPIR BH

3D cor LAVA ASPIR BH - 3D cor LAVA ASPIR BH

OTHERS

Protocol Notes



LAVA Considerations:
-Use LAVA for abdominal scanning, in particular, liver imaging
-LAVA automatically turns on the ARC Imaging Option. If you deselect ARC and select ASSET, you must have acquired a calibration scan prior to acquiring the LAVA scan. If any or all of these scan parameters are selected (ASSET, Hyperband, PURE), and if you select ON from the Calibration in Prescan menu, which is located on the Details tab, a calibration scan is acquired during Auto Prescan
Consider the following with Fat/Water SAT:
-The overall SNR of Improved Fat SAT images may slightly decrease.
Consider adjusting parameters with affect SNR, e.g. increase NEX, increase the FOV, decrease matrix size, decrease the Receive Bandwidth
-Uneven saturation can still occur as a result of local inhomogeneities, e.g. at air/tissue interfaces when the anatomy of interest is non-uniform
-Site or patient specific inhomogeneities may be unavoidable even at isocenter. The result can be uneven suppression.
Techniques for Fat Saturation include:
Chemical Fat SAT, STIR, ASPIR, and Flex.
-Consider using shom volumes to increase Fat SAT homogeneity.
Localized TG is used to improve the accuracy of the TG for a specifically designated area, such as cardiac, pelvis or prostate imaging, or when an anatomy is positioned at a patient's side. Localize TG can be used with any pulse sequence to improve signal homogeneity and is found under the Shim Volume menu

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

Ax LAVA Flex Pre - Ax LAVA Flex Pre	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	3D	Flip Angle	12
	Pulse Sequence	LAVA	TE	Min Full
	Imaging Options	EDR, Fast, Fluoro, ARC, Flex	Number of Echoes	2
	IDEAL	15	Receiver Bandwidth	62.50
	Phase	2.00	IMAGE ENHANCE	
	Slice	1.20	Filter Choice	A
	SCAN RANGE		USER CVS	
	FOV	42.0	User CV4	3.00
	Slice Thickness	2.30	User CV6	1.00
	Location per Slab	80	User CV30	1.00
	Overlap Locations	0	User CV34	70.00
	Number of Slices	1	User CV Mask2	8
	ACQ TIMING		MULTI-PHASE	
	Freq	288	Seperate Series	0
	Phase	200	Trigger Delay without AV	3
	Freq DIR	R/L	Mask Phase	0
	NEX	1.00	Mask Pause	0
	Phase FOV	0.80	Preserve	0
	Auto Shim	Auto	DIFFUSION	
	Phase Correction	No	Recon All Images	On
	RF Drive Mode	Single	# Synthetic b-values	1
	Excitation Mode	Selective	Synthetic b-value	1000.0;
	FMRI		CONTRAST	
	PSD Trigger	Internal	Contrast Yes/No	No
	View Order	Bottom/Up		
	# of Repetitions REST	0		
	# of Repetitions ACTIVE	0		
	SAT			
	Tag Type	None		
	TRICKS			
	Pause On/Off	On		
	Auto Subtract	0		
	Auto SCIC	2		

Ax LAVA Flex Pre - Ax LAVA Flex Pre

MSK-THR HIP NEW

Protocol: adult_abdomen_BOD -Gad Liver STANDARD aspir protocol - ARDL smartprep

Ax LAVA Flex Multiphase - Ax LAVA Flex Multiphase	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	3D	Flip Angle	12
	Pulse Sequence	LAVA	TE	Min Full
	Imaging Options	EDR, Fast, MPh, Fluoro, ARC, Flex	Number of Echoes	2
	IDEAL	3	Receiver Bandwidth	62.50
	Phase	2.00	IMAGE ENHANCE	
	Slice	1.20	Filter Choice	A
	SCAN RANGE		USER CVS	
	FOV	40.0	User CV4	3.00
	Slice Thickness	2.30	User CV6	1.00
	Location per Slab	80	User CV30	1.00
	Overlap Locations	0	User CV34	70.00
	Number of Slices	1	User CV Mask2	8
	ACQ TIMING		MULTI-PHASE	
	Freq	288	Slice per Location	3
	Phase	200	Delay after Acquisition	Minimum
	Freq DIR	R/L	Seperate Series	0
	NEX	1.00	Trigger Delay without AV	3
	# of Acq. Before Pause	1	Mask Phase	0
	Phase FOV	0.80	Mask Pause	0
	Auto Shim	Auto	DIFFUSION	
	Phase Correction	No	Recon All Images	On
	RF Drive Mode	Single	# Synthetic b-values	1
	Excitation Mode	Selective	Synthetic b-value	1000.0;
	FMRI		CONTRAST	
	PSD Trigger	Internal	Contrast Yes/No	Yes
	View Order	Bottom/Up		
	# of Repetitions REST	0		
	# of Repetitions ACTIVE	0		
	SAT			
	Tag Type	None		
	TRICKS			
	Pause On/Off	On		
	Auto Subtract	7		
	Auto SCIC	2		

Ax LAVA Flex Multiphase - Ax LAVA Flex Multiphase