

# MSK-THR HIP NEW

Protocol: adult\_lowerExtremities\_MSK -ARDL 2022-Knee std High Resolution

3-Plane Localizer Right - 3-Plane Localizer Right

IMAGING PARAMETERS	
Imaging Mode	2D
Pulse Sequence	Spin Echo
Imaging Options	Seq, EDR, TRF, Fast, SS, ARC
Phase	2.00
SCAN RANGE	
FOV	30.0
Slice Thickness	7.00
Slice Spacing	5.0
Overlap Locations	0
ACQ TIMING	
Freq	288
Phase	240
Freq DIR	Unswap
# of Acq. Before Pause	0
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective
FMRI	
PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off
OTHERS	
Protocol Notes	For accurate AIRx Knee calculations, ensure coverage spans above patella to below tibial plateau. Localizer must meet the following criteria: * Frequency and phase resolution >=192 or in-plane resolution of <=1.5 mm for both frequency and phase direction * Slice thickness <=10 mm * Slice coverage for all axial, sagittal and coronal plane >= 4 cm * Slice gap <=5 mm * TE for SSFSE <=80 ms

3-Plane Localizer Right - 3-Plane Localizer Right

SCAN TIMING	
TE	80.0
Number of Echoes	1
TR	600.0
Receiver Bandwidth	83.33
IMAGE ENHANCE	
Filter Choice	None
USER CVS	
User CV1	1.00
User CV39	1.00
User CV Mask2	256
MULTI-PHASE	
Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0
DIFFUSION	
Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;
CONTRAST	
Contrast Yes/No	No

3-Plane Localizer Right - 3-Plane Localizer Right

MSK-THR HIP NEW

Protocol: adult\_lowerExtremities\_MSK -ARDL 2022-Knee std High Resolution

SAG T2 FS - SAG T2 FS	<b>IMAGING PARAMETERS</b>		<b>SCAN TIMING</b>	
	Imaging Mode	2D	Flip Angle	150
	Pulse Sequence	FSE-XL	TE	68.0
	Imaging Options	EDR, TRF, Fast, ARC	Number of Echoes	1
	Phase	2.00	TR	5810.0
	<b>SCAN RANGE</b>		Echo Train Length	15
	FOV	16.0	Receiver Bandwidth	31.25
	Slice Thickness	3.00	<b>IMAGE ENHANCE</b>	
	Slice Spacing	0.5	Filter Choice	None
	Overlap Locations	0	<b>USER CVS</b>	
	Number of Slices	27	User CV12	1.00
	<b>ACQ TIMING</b>		User CV19	1.00
	Freq	384	User CV21	1.00
	Phase	320	TR Min	1500.0
	Freq DIR	A/P	TR Max	10000.0
	Fat Shift DIR	Normal (P)	User CV30	1.00
	NEX	1.00	User CV39	1.00
	Phase FOV	1.00	User CV Mask2	256
	Auto Shim	Auto	<b>MULTI-PHASE</b>	
	Phase Correction	Yes	Seperate Series	0
	RF Drive Mode	Single	Mask Phase	0
	Excitation Mode	Selective	Mask Pause	0
	<b>FMRI</b>		Preserve	0
	PSD Trigger	Internal	<b>DIFFUSION</b>	
	View Order	Bottom/Up	Recon All Images	On
	# of Repetitions REST	0	# Synthetic b-values	1
	# of Repetitions ACTIVE	0	Synthetic b-value	1000.0;
	<b>SAT</b>		<b>TRACKER</b>	
	SAT Location	S	Tracker Length	200.0
	SAT Location	I	Tracker Thickness	20.0
	Tag Type	None		
	Fat/Water Saturation	Fat Classic		
	<b>TRICKS</b>			
	Pause On/Off	On		
	Auto Subtract	0		
	Auto SCIC	3		
	<b>CONTRAST</b>			
	Contrast Yes/No	No		

SAG T2 FS - SAG T2 FS

## MSK-THR HIP NEW

Protocol: adult\_lowerExtremities\_MSK -ARDL 2022-Knee std High Resolution

**IMAGING PARAMETERS**

Imaging Mode	2D
Pulse Sequence	FSE-XL
Imaging Options	EDR, TRF, Fast, ARC
Phase	2.00

**SCAN RANGE**

FOV	16.0
Slice Thickness	3.00
Slice Spacing	0.5
Overlap Locations	0
Number of Slices	27

**ACQ TIMING**

Freq	384
Phase	320
Freq DIR	A/P
Fat Shift DIR	Normal (P)
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes
RF Drive Mode	Single
Excitation Mode	Selective

**FMRI**

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

**SAT**

SAT Location	S
SAT Location	I
Tag Type	None
Fat/Water Saturation	Fat Classic

**TRICKS**

Pause On/Off	On
Auto Subtract	0
Auto SCIC	3

**CONTRAST**

Contrast Yes/No	No
-----------------	----

**SCAN TIMING**

Flip Angle	150
TE	15.0
Number of Echoes	1
TR	1892.0
Echo Train Length	9
Receiver Bandwidth	31.25

**IMAGE ENHANCE**

Filter Choice	None
---------------	------

**USER CVS**

User CV12	1.00
User CV19	1.00
User CV21	1.00
TR Min	1500.0
TR Max	3500.0
User CV30	1.00
User CV39	1.00
User CV Mask2	256

**MULTI-PHASE**

Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0

**DIFFUSION**

Recon All Images	On
Multi b-values	1000.0;
Multi NEX Values	1.0;
# Synthetic b-values	1
Synthetic b-value	1000.0;

**TRACKER**

Tracker Length	200.0
Tracker Thickness	20.0

SAG PD FS - SAG PD FS

SAG PD FS - SAG PD FS

## MSK-THR HIP NEW

Protocol: adult\_lowerExtremities\_MSK -ARDL 2022-Knee std High Resolution

## IMAGING PARAMETERS

Imaging Mode	2D
Pulse Sequence	FSE-XL
Imaging Options	EDR, TRF, Fast, ARC
Phase	2.00

## SCAN RANGE

FOV	16.0
Slice Thickness	3.00
Slice Spacing	0.5
Overlap Locations	0
Number of Slices	34

## ACQ TIMING

Freq	384
Phase	320
Freq DIR	S/I
Fat Shift DIR	Normal (I)
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes
RF Drive Mode	Single
Excitation Mode	Selective

## FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

## SAT

Tag Type	None
Fat/Water Saturation	Fat Classic

## TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	3

## CONTRAST

Contrast Yes/No	No
-----------------	----

## SCAN TIMING

Flip Angle	150
TE	35.0
Number of Echoes	1
TR	2199.0
Echo Train Length	9
Receiver Bandwidth	31.25

## IMAGE ENHANCE

Filter Choice	None
---------------	------

## USER CVS

User CV12	1.00
User CV19	1.00
User CV21	1.00
TR Min	1500.0
TR Max	3500.0
User CV30	1.00
User CV39	1.00
User CV Mask2	256

## MULTI-PHASE

Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0

## DIFFUSION

Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;

## TRACKER

Tracker Length	200.0
Tracker Thickness	20.0

Cor PD FS - Cor PD FS

Cor PD FS - Cor PD FS

MSK-THR HIP NEW

Protocol: adult\_lowerExtremities\_MSK -ARDL 2022-Knee std High Resolution

Cor T1 - Cor T1

IMAGING PARAMETERS	
Imaging Mode	2D
Pulse Sequence	FSE-XL
Imaging Options	EDR, TRF, Fast, ARC
Phase	2.00
SCAN RANGE	
FOV	16.0
Slice Thickness	3.00
Slice Spacing	0.5
Overlap Locations	0
Number of Slices	22
ACQ TIMING	
Freq	384
Phase	352
Freq DIR	S/I
Fat Shift DIR	Normal (I)
NEX	2.00
# of Acq. Before Pause	0
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes
RF Drive Mode	Single
Excitation Mode	Selective
FMRI	
PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	3

SCAN TIMING	
Flip Angle	160
TE	Min Full
Number of Echoes	1
TR	471.0
Echo Train Length	4
Receiver Bandwidth	41.67
IMAGE ENHANCE	
Filter Choice	None
USER CVS	
User CV12	1.00
TR Min	400.0
TR Max	800.0
User CV30	1.00
User CV39	1.00
User CV Mask2	256
MULTI-PHASE	
Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0
DIFFUSION	
Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;
CONTRAST	
Contrast Yes/No	No

Cor T1 - Cor T1

## MSK-THR HIP NEW

Protocol: adult\_lowerExtremities\_MSK -ARDL 2022-Knee std High Resolution

## IMAGING PARAMETERS

Imaging Mode	2D
Pulse Sequence	FSE-XL
Imaging Options	EDR, TRF, Fast, ARC
Phase	2.00

## SCAN RANGE

FOV	16.0
Slice Thickness	3.00
Slice Spacing	0.5
Overlap Locations	0
Number of Slices	28

## ACQ TIMING

Freq	384
Phase	320
Freq DIR	A/P
Fat Shift DIR	Normal (A)
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes
RF Drive Mode	Single
Excitation Mode	Selective

## FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

## SAT

Tag Type	None
Fat/Water Saturation	Fat Classic

## TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	3

## CONTRAST

Contrast Yes/No	No
-----------------	----

## SCAN TIMING

Flip Angle	150
TE	34.0
Number of Echoes	1
TR	2012.0
Echo Train Length	9
Receiver Bandwidth	35.71

## IMAGE ENHANCE

Filter Choice	None
---------------	------

## USER CVS

User CV12	1.00
User CV19	1.00
User CV21	1.00
TR Min	1500.0
TR Max	3500.0
User CV30	1.00
User CV39	1.00
User CV Mask2	256

## MULTI-PHASE

Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0

## DIFFUSION

Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;

## TRACKER

Tracker Length	200.0
Tracker Thickness	20.0

Ax PD FS - Ax PD FS

Ax PD FS - Ax PD FS

## MSK-THR HIP NEW

Protocol: adult\_lowerExtremities\_MSK -ARDL 2022-Knee std High Resolution

Small FOV MR Bone - Small FOV MR Bone

IMAGING PARAMETERS	
Imaging Mode	3D
Pulse Sequence	SPGR
Imaging Options	Fast, SqP
PSD Name	ozteo
SCAN RANGE	
FOV	20.0
Slice Thickness	1.10
Overlap Locations	0
Number of Slices	1
ACQ TIMING	
Freq	200
Freq DIR	S/I
NEX	8.00
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Single
Excitation Mode	Selective
Spokes Per Segment	384
USER CVS	
User CV Mask2	0
MULTI-PHASE	
Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0
DIFFUSION	
Recon All Images	On
Multi b-values	1000.0;
Multi NEX Values	1.0;
# Synthetic b-values	1
Synthetic b-value	1000.0;
CONTRAST	
Contrast Yes/No	No

Small FOV MR Bone - Small FOV MR Bone

SCAN TIMING	
Flip Angle	2
Number of Echoes	1
Receiver Bandwidth	50.00
IMAGE ENHANCE	
Filter Choice	G
GATING/TRIGGER	
Pause After Navigator Prescan	0
FMRI	
PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	3

Small FOV MR Bone - Small FOV MR Bone

Small FOV MR Bone - Small FOV MR Bone